



# Maestro

SOFTWARE

## DESCRIPTION

Maestro is a software for managing patient data that can interface with all the Inventis device, controlling or getting exams results from them.

## OPERATING MODES

- Noah module
- Stand-Alone with proprietary database
- Stand-Alone without database functions (file based)
- Third-party Applications (exchange data with external applications): MII proprietary protocol and GDT standard

## AVAILABLE TESTS

- Audiology exams
- REM
- HIT
- Video-Otoscopy
- Video and Images acquisition
- Questionnaire

## COMPATIBLE DEVICES

- Triangle (Audiometer)
- Bell (Audiometer)
- Harp (Audiometer)
- Piano (Audiometer)
- Piccolo (Audiometer)
- Cello (Audiometer)
- Trumpet (REM, Audiometer)
- Drum (HIT)
- Flute (Middle-Ear Analyzer)
- Clarinet (Middle-Ear Analyzer)
- Viola (Middle-Ear Analyzer, Audiometer)
- Timpani (Handheld Tympanometer, Audiometer)
- Delfino (Otoscope)
- Harmonica (Otoscope)
- MediCam (Video-endoscopy camera)
- MediCam HD (HD Video-endoscopy camera)
- Flexia 2 (Video-Nasopharyngoscope)
- Analog-to-USB video converter
- Analog-to-USB video converter HD

## INTERFACE WITH EQUIPMENT

USB connection

## PC MINIMUM REQUIREMENTS

CPU: Intel® i5, 6<sup>th</sup> generation or above, 8 GB RAM  
 Operating System: Windows 10 and Windows 11, 64 bits  
 Display resolution: suggested 1920 x 1080 or higher (zoom factor 100% or 125%)  
 Local network (only for db installation in network folder): 1 Gb/s

## NOAH MINIMUM REQUIREMENTS

NOAH 4.4 or greater (or NOAH Engines supporting NOAH 4.4 or greater)

## SOFTWARE SUPPORT

Remote assistance tool: included

## SOFTWARE UPDATES

Available for free

## MII INTERFACE

Library for interfacing Maestro with 3<sup>rd</sup> party software  
 Available upon request (IT technician intervention requested)

## APPLICABLE STANDARDS

Software life cycle: IEC 62304  
 Health software: IEC 82304-1, IEC 81001-5-1

## CE CERTIFICATE

MDR 2017/745/EU Classification: Class IIa  
 Classification rule (Annex VIII, 2017/745): 11  
 Notified body: TÜV SÜD Product Service GmbH (0123)

Maestro is developed by:

### INVENTIS S.r.l.

CORSO STATI UNITI, 1/3  
 35127 PADOVA – ITALY  
 PHONE: +39.049.8962 844  
 FAX: +39.049.8966 343  
 info@inventis.it  
 www.inventis.it

Inventis® is a registered trademark of INVENTIS S.r.l.