



Maestro control panel

Accessory

DESCRIPTION

Maestro Control Panel (MCP) is an optional add-on to the Inventis audiology equipment, intended to ease the navigation and use of Maestro software, especially while performing audiometric exams.

INPUT/OUTPUT

Two set of connectors on the rear panel: one for cables coming from the audiometer/PC and the other one for connecting the detachable parts:

- Monitor headphone (3.5mm plug)
- Talk forward microphone (3.5mm plug)
- Desktop microphone (3.5mm plug)

OPERATOR INTERFACE

- 50 buttons (30 with built-in LEDs) with single specific audiometric function
- 4 configurable buttons
- 2 capacitive interrupters
- 2 knobs

COMPUTER INTERFACE

USB driverless connection

Compatible software: Inventis Maestro (from version 1.20.0)

POWER SUPPLY

USB powered

Power consumption: 5 Vdc, 0.25 A

MECHANICS

Size (WxDxH): 36 x 27 x 4 cm / 14.2 x 10.6 x 1.6 in

Weight: approx. 1.2 kg / 2.6 lbs

FREIGHT PACKING

Size (WxDxH): 45.5 x 38 x 8.5 cm / 17.9 x 15 x 3.3 in

Gross weight: approx. 2.5 kg / 5.5 lbs

APPLICABLE STANDARDS

Electrical safety: IEC 60601-1, Class II

EMC: IEC 60601-1-2

CE CERTIFICATE

MDR 2017/745/EU Classification: Class I

PRODUCT CODES

12080: Maestro control panel

INCLUDED PARTS

- USB connection cable
- 2 x stereo audio cables, 1.8m, 3.5mm plug - 3.5mm plug
- 1 x mono audio cable, 2m, 3.5mm plug - 6.35mm plug
- User manual

Maestro control panel is manufactured by:

INVENTIS S.r.l.

CORSO STATI UNITI, 1/3

35127 PADOVA – ITALY

PHONE: +39.049.8962 844

FAX: +39.049.8966 343

info@inventis.it

www.inventis.it

Follow us on LinkedIn

<https://www.linkedin.com/company/inventis-srl>

The Inventis Quality System complies with ISO 13485 standard.

Inventis® is a registered trademark of INVENTIS S.r.l.