



Maestro control panel

Accessory

DESCRIPTION

Maestro Control Panel (MCP) is an optional add-on to the Inventis audiology equipment, intended to ease the navigation and use of Maestro software, especially while performing audiometric exams. Even if compatible with most of the Inventis audiometric equipment, MCP is mostly recommended for being used with Cello and Trumpet audiometers.

RECOMMENDED FOR USE WITH

- Cello
- Trumpet AUD
- Trumpet REM & AUD
- Piccolo
- Harp, Piano (if controlled from the computer)

INPUT/OUTPUT

Two set of connectors on the rear panel: one for cables coming from the audiometer/PC and the other one for connecting the detachable parts:

- Monitor headphone (3.5mm plug)
- Talk forward microphone (3.5mm plug)
- Desktop microphone (3.5mm plug)

OPERATOR INTERFACE

- 50 buttons (30 with built-in LEDs) with single specific audiometric function
- 4 configurable buttons
- 2 capacitive interrupters
- 2 knobs

COMPUTER INTERFACE

USB driverless connection

Compatible software: Inventis Maestro (from version 1.20.0)

POWER SUPPLY

USB powered

Power consumption: 5 Vdc, 0.25 A

MECHANICS

Size (WxDxH): 36 x 27 x 4 cm / 14.2 x 10.6 x 1.6 in

Weight: approx. 1.2 kg / 2.6 lbs

FREIGHT PACKING

Size (WxDxH): 45.5 x 38 x 8.5 cm / 17.9 x 15 x 3.3 in

Gross weight: approx. 2.5 kg / 5.5 lbs

APPLICABLE STANDARDS

Electrical safety: IEC 60601-1, Class II

EMC: IEC 60601-1-2

CE CERTIFICATE

MDR 2017/745/EU Classification: Class I

PRODUCT CODES

12080: Maestro control panel

INCLUDED PARTS

- USB connection cable
- 2 x stereo audio cables, 1.8m, 3.5mm plug - 3.5mm plug
- 1 x mono audio cable, 2m, 3.5mm plug - 6.35mm plug
- Electronic IFU

Maestro control panel is manufactured by:

INVENTIS S.r.l.

VIA AUSTRIA, 38

35127 PADOVA – ITALY

PHONE: +39.049.8962 844

info@inventis.it

www.inventis.it

Follow us on LinkedIn

<https://www.linkedin.com/company/inventis-srl>

The Inventis Quality System complies with ISO 13485 standard.

Inventis® is a registered trademark of INVENTIS S.r.l.